

Steam Sterilizer XRD STE23 LCD

Advanced class B high technology steam sterilizer, the STE23 LCD has Dual Heating System (DHS), Multi Temperature Control, making the process more efficient than ever before.

Rapid USB port
Large Capacity water tank



→ Two types of record (Dual Recording)

1) USB

Master controller stores 2,000 fresh cycles; the data can be downloaded via USB. You can check and print all the records with a computer.

2) Micro-Printer(Optional) Prints cycle data when the cycle ends.




→ Big Capacity Open Tank

1) Open tank, easy-to-clean.


2) Built-in water quality sensor monitors water quality, reduces machine failure rate, and prolongs machine life time.




Short Sterilization Time



Constant Temperature Control



Heating Efficiency



Three-Dimensional Sterilization

Features

Rating Voltage	120 V
Power	1800 W
Capacity	23 L (6 gal)
Chamber Size	9.80" x 17.70"
External Dimensions	(WxHxD)18.18x18.11x23.62"
Net Weight	110 lbs