

trident



**(HD)** High Definition Digital Images  
I View captures images both horizontally and vertically, even in the most critical areas deep within the oral cavity, with a resolution of more than 20 visible line pairs per millimeter.

I-VIEW

xrddental

CMOS Intraoral Sensor

# I-VIEW

INTRAORAL DIGITAL SENSOR

Water-resistant  
Two ergonomically sizes  
Replaceable cable

I-VIEW is the CMOS intraoral sensor designed to capture highly detailed digital images in less time and great accuracy



#### DIRECT USB CONNECTION

Connect directly to a PC.  
Connection compatible with a standard 2.0 USB port.



#### LARGE ACTIVE SURFACE

1000 x 1500 pixel (20 x 30 mm) **size 1**  
1300 x 1700 pixel (26 x 34 mm) **size 2**



#### QUICK IMAGE CAPTURE

High quality images with the optimization of the job time



#### STRONG & DURABLE

Kevlar reinforced cable  
Replaceable cable




#### SAVE MONEY

We offer an affordable solution for the Dental Practice with high quality results

## THE BEST OF I-VIEW

Get instantaneously images whit less exposure to radiation.  
Increase technical and diagnostic capabilities.  
Organize and store thousands of images in your PC.



#### SOFTWARE DIRECT INTEGRATION


The I-View sensor natively works with DeepView, XrayVision and XV Lite and all the other brands under Apteryx including: Belmont XV, Cliniview, Dentalind, Oriview, Dentimax advanced, XL Dent, Image XL, LynxVision Pro, MaxiViewer, Benco OneView, Prof. Suni, QSI Image, Sigma Image XRV.



#### TWAIN COMPATIBLE

Works in all leading imaging software available in the market: Dexis V10 and newer, TigerView, eaglesoft Advanced Imaging, Dentrix Ascend, Curve Dental, Mediadent, Patient Gallery, Dovetail, cloud9software, TDO, Cleardent and more.

**Note:** The Software and Twain Driver are not included with the product. We recommend Apteryx brand.



# xrddental

Dental Equipment - [www.xrddental.com](http://www.xrddental.com)

XRD Dental Equipment - Orlando FL 32836  
+ 1 (407) 319-1348 - [info@xrddental.com](mailto:info@xrddental.com)