

# I-View Gold VET

## CMOS Intraoral Sensor

Practical cover to protect the active area from scratches and damages.

Available in Sizes #2 and # 1

Water-resistant



## Two Sizes

Thinner body and even surface  
Rounded corners for a comfortable and easy positioning



### Dimensions

Size 1.0 25.4 x 36.8mm

Size 2.0 31.3 x 42.9mm



### Specifications

Detector CMOS

Pixel size 18 Micron


Theoretical Resolution 28 lp/mm

Dynamic Range 12 bit



### Deep View 4® Software

New generation of Deep-View software, designed to efficiently create, organize store and share, your pets' data files and images.



## Direct USB Connection



- Easy images acquisition device configuration
- Advanced image filters
- Wide choice and customization of templates for Full Mouth Series exams
- Voice guidance option for image acquisition
- Images import/export in DICOM and other graphic formats
- Printing images with layout customization
- TWAIN included for integration with other softwares
- Pre-Calibrated and easy installation. Plug & Use



# xrddental

Dental Equipment - [www.xrddental.com](http://www.xrddental.com)

XRD Dental Equipment - Orlando FL 32836  
+ 1 (407) 319-1348 - [info@xrddental.com](mailto:info@xrddental.com)